

# ULTRASONIC BATH / CLEANER

International standard ultrasonic bath / cleaning machine **88** series



## Features:

- The most popular design
- This Ideal equipment for laboratory, operating room, endoscope room, and sterilization room
- U.S. control system operation convenience
- Independent degassing function

## Technical Specifications:

Model	External Dimensions L x W x H (mm)	Internal Dimensions L x W x H (mm)	Capacity (L)	Frequency (kHz)	Power (W)	Cover	Body	Timing/Heating	Drainage	Basket
KS-8892	406 x 305 x 368	300 x 150 x 150	6	40	200	Yes	Plastic	Digital display/yes	Yes	Yes
KS-8893	406 x 394 x 368	300 x 240 x 150	10	40	200	Yes	Plastic	Digital display/yes	Yes	Yes
KS-8895	606 x 455 x 368	500 x 300 x 150	22	40	500	Yes	Plastic	Digital display/yes	Yes	Yes

International standard ultrasonic bath cleaning machine 99 series



**Features:**

- Low noise
- Equipped with net basket and pipe
- U.S. control system operation convenience
- The exterior is completely designed in accordance to export standards

**Single Frequency Technical Specifications:**

Model	Power (W)	Frequency (kHz)	External Dimensions L x W x H (mm)	Internal Dimensions L x W x H (mm)	Capacity (L)	Drainage	Basket	Timing/Heating	Degassing Function
9980	80	40	230 x 210 x 245	150 x 140 x 100	2	No	Yes	Digital Display	Yes
9981	120	40	310 x 210 x 245	240 x 140 x 100	3	No	Yes	Digital Display	Yes
9982	180	40	370 x 225 x 300	300 x 150 x 150	6	Yes	Yes	Digital Display	Yes
9983	300	40	370 x 310 x 300	300 x 240 x 150	10	Yes	Yes	Digital Display	Yes
9985	600	40	590 x 380 x 300	500 x 300 x 150	22	Yes	Yes	Digital Display	Yes
9986	800	40	590 x 380 x 350	500 x 300 x 200	30	Yes	Yes	Digital Display	Yes

**Double Frequency Technical Specifications:**

Model	Power (W)	Frequency (kHz)	External Dimensions L x W x H (mm)	Internal Dimensions L x W x H (mm)	Capacity (L)	Drainage	Basket	Timing/Heating
9980-2	80	40/59 Double Frequency	230 x 210 x 245	150 x 140 x 100	2	No	Yes	Digital Display
9981-2	120	40/59 Double Frequency	310 x 210 x 245	240 x 140 x 100	3	No	Yes	Digital Display
9982-2	180	40/59 Double Frequency	370 x 225 x 300	300 x 150 x 150	6	Yes	Yes	Digital Display
9983-2	300	40/59 Double Frequency	370 x 310 x 300	300 x 240 x 150	10	Yes	Yes	Digital Display
9985-2	600	40/59 Double Frequency	590 x 380 x 300	500 x 300 x 150	22	Yes	Yes	Digital Display
9986-2	800	40/59 Double Frequency	590 x 380 x 350	500 x 300 x 200	30	Yes	Yes	Digital Display

## Table ultrasonic bath / cleaning machine **EI/EII** series



### Features:

- Low noise
- Equipped with net basket and pipe
- U.S. control system operation convenience
- The exterior is completely designed in accordance to export standards
- Independent degassing function

### EI Technical Specifications:

Model	External Dimensions L x W x H (mm)	Internal Dimensions L x W x H (mm)	Capacity (L)	Frequency (kHz)	Power (W)	Timing/ Heating	Drainage	Cover	Body	Basket
KS-120EI	260 x 160 x 240	240 x 140 x 100	3	40	120	Digital display	No	Yes	Stainless Steel	Yes
KS-150EI	320 x 175 x 240	300 x 150 x 100	4	40	150	Digital display	No	Yes	Stainless Steel	Yes
KS-180EI	350 x 175 x 310	300 x 150 x 150	6	40	180	Digital display	Yes	Yes	Stainless Steel	Yes
KS-250EI	350 x 264 x 310	300 x 240 x 150	10	40	250	Digital display	Yes	Yes	Stainless Steel	Yes
KS-300EI	350 x 264 x 310	300 x 240 x 150	10	40	300	Digital display	Yes	Yes	Stainless Steel	Yes
KS-500EI	550 x 320 x 310	500 x 300 x 150	22	40	500	Digital display	Yes	Yes	Stainless Steel	Yes
KS-600EI	550 x 320 x 310	500 x 300 x 150	22	40	600	Digital display	Yes	Yes	Stainless Steel	Yes

### EII Technical Specifications:

Model	External Dimensions L x W x H (mm)	Internal Dimensions L x W x H (mm)	Capacity (L)	Frequency (kHz)	Power (W)	Timing/ Heating	Drainage	Cover	Body	Basket
KS-120EII	260 x 160 x 240	240 x 140 x 100	3	40/59	120	Digital display	No	Yes	Stainless Steel	Yes
KS-150EII	320 x 175 x 240	300 x 150 x 100	4	40/59	150	Digital display	No	Yes	Stainless Steel	Yes
KS-180EII	350 x 175 x 310	300 x 150 x 150	6	40/59	180	Digital display	Yes	Yes	Stainless Steel	Yes
KS-250EII	350 x 264 x 310	300 x 240 x 150	10	40/59	250	Digital display	Yes	Yes	Stainless Steel	Yes
KS-300EII	350 x 264 x 310	300 x 240 x 150	10	40/59	300	Digital display	Yes	Yes	Stainless Steel	Yes
KS-500EII	550 x 320 x 310	500 x 300 x 150	22	40/59 (28/40)	500	Digital display	Yes	Yes	Stainless Steel	Yes
KS-600EII	550 x 320 x 310	500 x 300 x 150	22	40/59 (28/40)	600	Digital display	Yes	Yes	Stainless Steel	Yes

**NovaScientific**

**NovaScientific Resources (M) Sdn. Bhd.**

No. 12A-2A, Jalan PJU 1/3B, Sunwaymas Commercial Centre, 47301 Petaling Jaya,  
Selangor Darul Ehsan, Malaysia

Tel: 03-7805 5766 Fax: 03-7805 5866

E-mail: [novascientific@gmail.com](mailto:novascientific@gmail.com), Website: [www.novascientific.com.my](http://www.novascientific.com.my)